



# Rotaid AED Cabinets

## Accessoires - Totem



English  
version

# Product sheet

*The Rotaid AED Cabinets totem allows for a stand alone in- and outdoor installation of your defibrillator. The solid construction and high quality appearance make sure your defibrillator stands out and can be securely placed in practically all environments and situations.*

The modular set-up ensures an easy installation and the totem can easily be demounted when used for training and trade show purposes.

On the front, back and top of the totem, information and branding panels can be mounted. These panels are custom made, connecting to the customer's request and location, making sure the totem aligns perfectly with its surroundings.

The cabinet is safely mounted on the totem and the solid steel construction allows for installation in the most demanding environments. The double powder coated finish ensures a long-lasting appearance.

## Appliances

- **Training and fares**
- **Stand alone installations**
- **Locations with no wall mounting possibilities**
- **Outdoor installations**
- **Renting and eventing**

## Design

- **Custom made design**
- **Full colour UV protected printing**
- **Allows for fixed floor mounting**
- **In combination with all Rotaid cabinets**

- ✦ For indoor and outdoor usage
- ✦ Solid and safe construction
- ✦ Storage box for transformer unit
- ✦ Custom made designs and panels



# Specifications

## Dimensions

### Total

Width: 700 mm  
Height: 1615 mm  
Weight: 40 kg

### Footplate

Width: 700 mm round  
Thickness: 5 mm high  
Material: Steel, powder-coated

### Uprights / Stairs

Width: 560 mm  
Height: 1715 mm  
Thickness pole: 60 mm  
Material: Steel, powder-coated

### Panel

Width: 394 mm  
Height: 1075 mm  
Material: Dibond, 3 mm thick

## Article numbers

### 30.9000

Rotaid Totem with foot and front sign.

### 30.9000 + 30.0001

Rotaid Totem with foot, front sign and backboard.

### 30.9001

Rotaid Totem with foot, front sign and sunroof.

## Modules

### 1. Sunroof

The perfect protection to maximize shade.

### 2. Backboard (Custom made)

Provide additional guidance and instructions.

### 3. Panel (Custom made)

Promote your totem with a custom made sign.

### 4. Additional color

Available in every color.



**AED**

AUTOMATISCHE EXTERNE DEFIBRILLATOR



**ROTAID™**  
AED CABINETS



**Article number - 30.9000**  
Rotaid Totem with foot and front sign.



**Article number - 30.9000+30.0001**  
Rotaid Totem with foot, front sign and backboard.







**Article number - 30.9001**

Rotaid Totem with foot, front sign and sunroof.



---

 [rotaid.com](http://rotaid.com)    [info@rotaid.com](mailto:info@rotaid.com)    +31 (0)478 769 052  
 Keizersveld 52D   5803 AN, Venray   The Netherlands